

# Zirlux™ Shading Liquids

Zirlux Shading Liquids create highly esthetic, natural looking restorations that blend seamlessly into the surrounding dentition. Paired with the Zirlux 16+ White zirconia, these liquids create accurate and consistent results that you can rely on to produce all 16 VITA®\* shades.

Characterize your monolithic units for more detailed customization using the Zirlux FC2 Characterization stains post sintering.

### Benefits:

- Available in 16 VITA® dentin shades, 5 incisal, and 11 effect shades
- Apply to restorations in the green state and sinter according to manufacturer's instructions for use
- Can be used with any post staining system designed for zirconia



### Dentine Kit

**(901-2419)** 30mL **\$399.99**

Contains all 16 Dentine shades.



Dentine	Shade
901-2374	OM1
901-2375	OM2
901-2377	OM3
901-2380	A1
901-2381	A2
901-2382	A3
901-2383	A3.5
901-2384	A4
901-2385	B1
901-2386	B2
901-2387	B3
901-2388	B4
901-2389	C1
901-2390	C2
901-2391	C3
901-2392	C4
901-2394	D2
901-2395	D3
901-2396	D4
<b>100mL</b>	<b>\$58.99</b>



Effects	Shade
901-2403	Light Pink
901-2405	Dark Pink
901-2406	Purple
901-2408	Yellow
901-2409	Orange
901-2410	Light Brown
901-2411	Dark Brown
901-2412	Blue
901-2413	Light Grey
901-2414	Dark Grey
901-2415	Graphite
<b>30mL</b>	<b>\$57.99</b>



Incisal	Shade
901-2397	1
901-2398	2
901-2400	3
901-2401	4
901-2402	5
<b>15mL</b>	<b>\$13.99</b>



Dentin	Shade
901-2416	Dentin Lightener
<b>30mL</b>	<b>\$58.99</b>

### Effects Shades



\*VITA is a registered trademark of Vita Zahnfabrik.