

Zirlux[®] Temp Multi

The esthetic solution
for *temporary restorations*.

Dental laboratories can use Zirlux Temp Multi to fabricate temporary dental restorations for a wide range of cases, including more complex indications. Zirlux Temp Multi has a gradient transition of layers throughout the disc making it an esthetic solution with high color stability which easily polishes to a high luster.



Benefits:

- Available in 5 shades—Bleach, A1, A2, A3, B1 and B2
- Compatible with industry-leading milling machines

Indications:

- Temporary Crowns
- Temporary Bridges

Material Properties

Flexural Strength: >100MPa
Solubility: <0.8 µg/mm³
Water uptake: <20 µg/mm³
Modulus of elasticity: >2,200 MPa



Shade	98.5mm		89 x 71mm		95mm		100mm	
A1	98 x 12mm	900-7848	89 x 71 x 13mm	900-7868	95 x 16mm	900-7860	100 x 16mm	900-7876
	98 x 16mm	900-7849	89 x 71 x 20mm	900-7869	95 x 20mm	900-7861	100 x 20mm	900-7877
	98 x 20mm	900-7850	—	—	—	—	—	—
A2	98 x 12mm	900-7851	89 x 71 x 13mm	900-7870	95 x 16mm	900-7862	100 x 16mm	900-7878
	98 x 16mm	900-7852	89 x 71 x 20mm	900-7871	95 x 20mm	900-7863	100 x 20mm	900-7879
	98 x 20mm	900-7853	—	—	—	—	—	—
A3	98 x 12mm	900-7854	89 x 71 x 13mm	900-7872	95 x 16mm	900-7864	100 x 16mm	900-7880
	98 x 16mm	900-7855	89 x 71 x 20mm	900-7873	95 x 20mm	900-7865	100 x 20mm	900-7881
	98 x 20mm	900-7856	—	—	—	—	—	—
B2	98 x 12mm	900-7857	89 x 71 x 13mm	900-7874	95 x 16mm	900-7866	100 x 16mm	900-7882
	98 x 16mm	900-7858	89 x 71 x 20mm	900-7875	95 x 20mm	900-7867	100 x 20mm	900-7883
	98 x 20mm	900-7859	—	—	—	—	—	—
BL	98 x 12mm	900-7845	—	—	—	—	—	—
	98 x 16mm	900-7846	—	—	—	—	—	—
	98 x 20mm	900-7847	—	—	—	—	—	—
		\$125.99		\$125.99		\$125.99		\$125.99