

# Zirlux<sup>®</sup> Temp



Zirlux Temp is a monochromatic PMMA material used to fabricate temporary restorations. It is available in all 16 Vita shades and a variety of sizes for compatibility with leading milling machines. Zirlux Temp can be used to check the fit and function of a restoration on a model with ease.

#### Benefits:

- Easy to polish to a high luster
- High color stability
- Can be used with or without a veneering composite to achieve higher esthetics
- Compatible with all temporary cements

#### Indications:

- Temporary Crowns
- Temporary Bridges

#### Material Properties

Flexural Strength: >100MPa  
 Solubility: <0.8 µg/mm<sup>3</sup>  
 Water uptake: <20 µg/mm<sup>3</sup>  
 Modulus of elasticity: >2,200 MPa



Shade	98.5 x 18mm	95 x 20mm	71 x 20mm	100 x 20mm
A0	900-7899	—	—	—
A1	900-7900	900-7916	900-7920	900-7924
A2	900-7901	900-7917	900-7921	900-7925
A3	900-7902	900-7918	900-7922	900-7926
A3.5	900-7903	—	—	—
A4	900-7904	—	—	—
B1	900-7905	—	—	—
B2	900-7906	900-7919	900-7923	900-7927
B3	900-7907	—	—	—
B4	900-7908	—	—	—
C1	900-7909	—	—	—
C2	900-7910	—	—	—
C3	900-7911	—	—	—
C4	900-7912	—	—	—
D3	900-7914	—	—	—
D4	900-7915	—	—	—
	\$68.99	\$68.99	\$68.99	\$68.99