




Zirlux® 3D Print Resins

Zirlux® 3D Print Resins are designed to fit seamlessly into your digital workflow. We've validated these resins on leading 3D printers to ensure consistent, high-quality results.

Use the chart below to find your printer and select the appropriate print profile.

- Look for the  to identify official Zirlux print profiles.
- For printers where a Zirlux® profile is not yet available, select the KeyPrint® print profile with confidence to produce the same reliable results.



-  Zirlux Print Profiles available
-  KeyPrint Print Profile

	Zirlux® Model Print Resin	Zirlux® Gingiva Mask Print Resin	Zirlux® Surgical Guides Print Resin	Zirlux® Nightguard & Splint Print Resin - Flexible	Zirlux® Nightguard & Splint Print Resin - Rigid
Asiga® MAX (385nm)					
Asiga® MAX 2 (385nm)					
Asiga® Ultra 50 (385nm)					
Asiga® PRO 4k (385nm)					
Ackuretta FreeShape (405nm)					
Ackuretta Dentiq (405nm)					
Ackuretta SOL (405nm)					
Ackuretta SOL Plus (405nm)					
Anycubic Photon Mono X 6Ks (405nm)					
Anycubic Photon Mono M5 (405nm)					
B9 Creation Dent XL (385nm)					
Carbon M2 (385nm)					
Carbon M3 (385nm)					
Carbon M3 Max (385nm)					
Dekema trix™ Print (385nm)					
Elegoo Mars 3 Pro 4k (385nm)					
EPAX DX10 Pro (405nm)					
Flashforge Hunter (405nm)					
Formlabs Form 4B					
HeyGears A2D (385nm)					
HeyGears ChairSide Pro (385nm)					
HeyGears Ultracraft A3D (385nm)					
Microlay Eve (405nm)					
Microlay Versus (385nm)					
MiiCraft 125 (405nm)					
MiiCraft Ultra (405nm)					
Nexa3D NXD200 (405nm)					
Nexa3D XiP (405nm)					
Phrozen Sonic XL 4k (405nm)					
Phrozen Sonic Mini 4k (405nm)					
Planmeca Creo C5 (405nm)					
Rapid Shape D20+, D30+, D20II, D30II, D40II (385nm)					
Rapid Shape D50+, D70, D90, D100 (385nm)					
RayShape® Shape 1 (405nm)					
Roland® 3DXPRINT					
SHINING 3D® AccuFab-D1/D1s (405nm)					
SHINING 3D® Accufab-L4D (405nm)					
Sisma Everes UNO (405nm)					
Stratasys Origin One Dental (385nm)					
Structo ST-01					
W2P SolFlex 135/150/163/170/350/363/650 (385nm)					
W2P SolLab 335/370/450 (385nm)					
Zortrax Inkspire 2 Clinical (405nm)					
ZyloDent					

Specially formulated for Zirlux®, tested and validated using KeyPrint® settings and workflows for seamless integration.



Ready to optimize your workflow?

Visit www.zirlux.com/printresins or scan QR code to learn more.

March 2026