

Zirlux® 3D Print Resins

Zirlux® 3D Print Resins are designed to fit seamlessly into your digital workflow. We’ve validated these resins on leading 3D printers to ensure consistent, high-quality results.

Use the chart below to find your printer and select the appropriate print profile.

- Look for the  to identify official Zirlux print profiles.
- For printers where a Zirlux® profile is not yet available, select the KeyPrint® print profile with confidence to produce the same reliable results.

					
	Zirlux® Model Print Resin	Zirlux® Gingiva Mask Print Resin	Zirlux® Surgical Guides Print Resin	Zirlux® Nightguard & Splint Print Resin - Flexible	Zirlux® Nightguard & Splint Print Resin - Rigid
Asiga® MAX (385nm)					
Asiga® Ultra 50 (385nm)					
Asiga® PRO 4k (385nm)					
Ackuretta FreeShape (405nm)					
Ackuretta Dentiq (405nm)					
Ackuretta SOL (405nm)					
Ackuretta SOL Plus (405nm)					
Anycubic Photon Mono X 6Ks (405nm)	✓			✓	✓
Anycubic Photon Mono M5 (405nm)	✓		✓	✓	✓
B9 Creation Dent XL (385nm)		✓	✓		
Carbon M2 (385nm)					
Carbon M3 (385nm)					
Carbon M3 Max (385nm)					
Dekema trix™ Print (385nm)	✓	✓	✓		
Elegoo Mars 3 Pro 4k (385nm)	✓	✓	✓	✓	✓
EPAX DX10 Pro (405nm)			✓		
Flashforge Hunter (405nm)		✓	✓		
HeyGears A2D (385nm)				✓	✓
HeyGears ChairSide Pro (385nm)					✓
HeyGears Ultracraft A3D (385nm)				✓	✓
Microlay Eve (405nm)	✓	✓	✓	✓	
Microlay Versus (385nm)	✓	✓	✓	✓	
MiiCraft 125 (405nm)	✓	✓	✓		
MiiCraft Ultra (405nm)	✓	✓	✓		
Nexa3D NXD200 (405nm)	✓		✓	✓	✓
Nexa3D XiP (405nm)	✓	✓	✓	✓	
Phrozen Sonic XL 4k (405nm)	✓	✓	✓	✓	✓
Phrozen Sonic Mini 4k (405nm)	✓	✓	✓	✓	✓
Planmeca Creo C5 (405nm)	✓		✓	✓	
Rapid Shape D20+, D30+, D20II, D30II, D40II (385nm)					
Rapid Shape D50+, D70, D90, D100 (385nm)	✓				
RayShape® Shape 1 (405nm)	✓		✓	✓	✓
Roland® 3DXPRINT	✓		✓	✓	✓
SHINING 3D® AccuFab-D1/D1s (405nm)					
SHINING 3D® Accufab-L4D (405nm)					
Sisma Everes UNO (405nm)	✓	✓	✓	✓	✓
SprintRay Pro (405nm)		✓	✓	✓	✓
SprintRay Pro 2 (385nm)				✓	
MoonRay® S100/D75 (385nm)			✓		
Stratasys Origin One Dental (385nm)		✓	✓	✓	✓
Structo ST-01	✓				
W2P SolFlex 135/150/163170/350/363/650 (385nm)		✓	✓		
W2P SolLab 335/370/450 (385nm)		✓	✓		
Zortrax Inkspire 2 Clinical (405nm)	✓		✓	✓	

Specially formulated for Zirlux®, tested and validated using KeyPrint® settings and workflows for seamless integration.



Ready to optimize your workflow?
Visit www.zirlux.com/printresins or scan QR code to learn more.

October 2025